

Connect Your Scanner with Sun Dental Labs UK

Free Shipping, Free Implant Scan Bodies

Free digital Models

3Shape

- Visit <https://us.3shapecommunicate.com>
- Select "CONNECTIONS"
- Request a connection with Sun Dental Labs UK–email address eucad@sundentallabs.com
- After your connection is approved, select Sun Dental Labs when sending files

iTero

- Please contact your iTero representative and request a link to Sun Dental Labs UK– you can use email address info@sundentallabs.co.uk

Dentsply Sirona / DS Core

- Log in to your account using the 'Connect Case Centre' or Ds Core software
- Choose 'My Favourite Recipients'
- Search for Sun Dental Labs UK and click add – you can use email address eucad@sundentallabs.com

Carestream / Dexis IS connect

- Visit "Carestream " or 'Dexis IS connect' on your scanner
- Search for Sun Dental Labs UK or you can use email address to find us - cadcam@sundentallabs.co.uk
- Add Sun Dental Labs UK as a production lab
- Select Sun Dental Labs UK when submitting scans

Shinning 3D

- Visit Shining 3D online portal within the "Connection" tab – add Sun Dental Labs UK as a production lab – or find us using email address info@sundentallabs.co.uk

Medit

- Log into www.meditlink.com with the administrator account
- Go to 'Partners '
- Search for Sun Dental Labs UK or find us using email address info@sundentallabs.co.uk and click the 'Request Partnership' button
- You'll see that the Sun Dental Labs UK has been added to the "My Partner List"

For other scanners please send files to info@sundentallabs.co.uk via wetransfer.com or get in touch with us and we will install and create an account for the scanner that you are using .

Do you have a question about our products or services? Need help on submitting a case or creating a MySun Account? Our team is happy to help! Simply give us a call to talk to a friendly Sun Dental Customer Service Representative on 01226 786163 or visit our website

<http://www.sundentallabs.co.uk/>

The Future of Dentistry is Here

We believe it is time for a better patient experience. We believe now is the time for dentistry to go digital.

